


Shape Your Waist

Shaken Massager

PEMF Smart Fat-Burning Device





**Do you want to have
a nice body shape ?**

Measure Your Waist



Waist Measurement Method:

1. Use a soft tape measure with 0.1cm precision.
2. Wear lightweight clothing, stand upright, look forward, and relax your abdomen.
3. Wrap the tape measure around the thinnest part of your waist (1-2cm above the navel) or directly around the navel.
4. Ensure the tape measure is snug against the skin, with enough room to fit two fingers comfortably.
5. Record the measurement to the nearest 0.1cm.

World Health Organization (WHO) Standards for Abdominal Obesity:


International Standard (applicable to Europe, Africa, and other regions):

Category	Male	Female
Normal	<94cm	<80cm
Abdominal Obesity	≥94cm	≥80cm
High Risk	≥102cm	≥88cm


Asian Population Standard (applicable to East Asia, Southeast Asia)

Category	Male	Female
Normal	<90cm	<80cm
Abdominal Obesity	≥90cm	≥80cm
High Risk	≥95cm	≥90cm

Abdominal Obesity



Peripheral Obesity



What are the health risks of abdominal obesity?


High blood pressure

High blood lipids

Type 2 diabetes

Non-alcoholic fatty liver disease

Osteoarthritis



Stroke

Sleep apnea

Coronary heart disease

Heart failure

increased cancer risk



**Common weight loss struggles
today:
too busy, too lazy, cravings too
strong, no results**



To address
the uncertain
distribution of fat
through
established physical laws

Shaken Massager PEMF Intelligent Fat Removal Device



Smart Technology Leaves Fat Nowhere to Hide

Shaken Massager
PEMF Smart Fat-Burning Device



Ultrasou

Fat Reduction

Vibration

Body Shaping

Benefit belly internal organs

7.PEMF


Red Light

Heat Therapy

Massage

RF (RadioFrequency)

Fat Reduction: Intelligent Sonic Technology to Target Stubborn Fat




Ultrasonic

Utilizing CellBreak™ Cell Disruption Technology, it precisely targets the subcutaneous lipid layer (10-20mm deep), employs high frequency resonance (1MHz) to precisely induce resonant dissociation of adipocyte membranes, and catalyzes the breakdown of lipids into



Vibration



High-frequency (TBD) vibrations stimulate the contraction of waist and abdominal muscles, compressing fat cells.

Heat Therapy



Thermal penetration awakens abdominal tissues and cells, promotes blood circulation, and accelerates fat metabolism.

Focused ultrasound destroys deep-layer fat cells.

Cavitation Effect


The pressure changes created by the ultrasound waves form tiny bubbles in the water. These bubbles repeatedly expand and contract with the sound wave vibrations. When the sound wave energy reaches a certain intensity, the bubbles collapse violently, releasing powerful shockwaves and microjets. This directly tears the fat cell membranes, causing triglycerides to leak into the intercellular space.

Mechanical Shearing Action

The high-frequency vibrations of ultrasound generate continuous mechanical shear forces on the fat tissue. These forces act directly on the fat cell membranes. Under repeated mechanical stretching, the membranes experience increased permeability and structural damage, ultimately leading to their rupture.

Focused ultrasound destroys deep-layer fat cells


Before



After



Subject's
Waist
Circumference



—From "A Study on the Efficacy and Safety of Focused Ultrasound for Fat Reduction and Body Contouring," Peking Union Medical College Hospital



Original Article

Leptin level decreases after treatment with the combination of Radiofrequency and Ultrasound cavitation in response to the reduction in adiposity

Mahla Arabpour-Dahoue ^{a,b,1}, Elham Mohammadzadeh ^{f,1}, Amir Avan ^f, Pouya Nezafati ^c, Samira Nasrfard ^b, Hamideh Ghazizadeh ^f, Mehrane Mehramiz ^b, Mohammad Safarian ^a, Mohsen Nematy ^a, Lida Jarahi ^d, Gordon A. Ferns ^e, Abdolreza Norouzy ^{a,*}, Majid Ghayour-Mobarhan ^{f,**}

^a Biochemistry of Nutrition Research Center, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran

^b Department of Nutrition, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran

^c Student Research Committee, Mashhad University of Medical Sciences, Mashhad, Iran

^d Department of Community Medicine, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran

^e Brighton & Sussex Medical School, Division of Medical Education, Falmer, Brighton, Sussex, BN1 9PH, UK

^f Metabolic Syndrome Research Center, Mashhad University of Medical Sciences, Mashhad, Iran

ARTICLE INFO

Article history:

Received 11 January 2019

Accepted 22 January 2019

Keywords:

Leptin

Obesity (overweight)

Radiofrequency

Ultrasound

References

- «Leptin level decreases after treatment with the combination of Radiofrequency and Ultrasound cavitation in response to the reduction in adiposity»
- «Treatment by Combination of Non-Invasive Radio-Frequency and Ultrasound Cavitation Ameliorates the Condition of Obesity»
- «Impact of Contactless Apoptosis-Inducing RF on Temperature of Human Skin Surface and Subcutaneous Layer as well as Porcine Histology: A Pilot Study»
- «Radiofrequency Heating and High-Intensity Focused Electromagnetic Treatment Delivered Simultaneously: The First Sham-Controlled Randomized Trial»
- «Induction of Fat Apoptosis by a Non-Thermal Device: Mechanism of Action of Non-Invasive High-Intensity Electromagnetic Technology in a Porcine Model»

Mod Care J, 2019 July; 16(3):e90769.

Published online 2019 June 16.



Treatment by Combination of Non-Invasive Radio-Frequency and Ultrasound Cavitation Ameliorates the Condition of Obesity

Mina Safari Bidokhti¹, Zahra Khorasanchi¹, Elham Moh Naseri³, Hamideh Ghazizadeh², Maryam Shahi², Mehr Afshari¹, Hamideh Safarian², Amir Avan² and Majid Gh

¹Department of Nutrition, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran
²Metabolic Syndrome Research Center, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran
³Student Research Committee, Mashhad University of Medical Sciences, Mashhad, Iran
⁴Cardiac Rehabilitation Research Center, Cardiovascular Research Institute, Mahan University

*Corresponding author: Metabolic Syndrome Research Center, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran

Received 2019 February 17; Revised 2019 May 17; Accepted 2019 June 03.

Abstract

Background: Different non-invasive body contouring techniques are radiofrequency (RF), and ultrasound cavitation (USC).

Objectives: The aim of the present study was to evaluate the effects of combination of RF and USC on body contouring.

Methods: One hundred twenty-nine obese and overweight participants (group 1), abdomen and hips (group 2), abdomen and thighs (group 3), and thighs (group 4). The participants received RF/USC twice a week for five weeks.

Results: In the results, except for hips circumferences we observed a significant decrease in all other measurements after treatment with combined RF and USC therapy. Although, there was no significant difference in body weight before and after treatment.

Conclusions: Our findings demonstrated that combined use of RF and USC can reduce abdomen, waist, and thigh circumferences as well as improve body contouring.

Keywords: Radiofrequency, Ultrasound Cavitation, Obesity/Overweight

ABSTRACT

Background: Obesity and overweight are major public health problem. Different strategies have been developed for body contouring including Radiofrequency (RF) and Ultrasound (US). The aim of this study was to investigate changes in serum-leptin as a potential-modulator of food/energy intake, in overweight-women receiving RF/US and diet-therapy as well as the effect of therapy on modulation of lipid-profile and body-fat-mass.

Methods: Fifty overweight-females were enrolled in current randomized-clinical-trial and randomly divided into two groups. The case group received RF/US twice a week for 5 weeks with a low calorie diet whilst the control-group received just a low calorie diet. Demographic, biochemical markers as well as anthropometric measurements were determined.

Results: The level of leptin was reduced from 1.29 ± 0.32 ng/ml to 1.14 ± 0.34 ng/ml in case group, before and after treatment. No significant difference was observed in the serum-leptin levels of subjects in the control group. The combination of RF and US decreased the leptin-level. In particular, the combination of RF and US decreased the leptin-level significantly more than the control group.

Conclusions: Our findings demonstrated the reduction of the leptin after treatment with the combination of RF and US. This reduction was inversely correlated with LDL levels. Weight was reduced in case and control groups. The combination of RF and US decreased the leptin-level significantly more than the control group.

Keywords: Radiofrequency, Ultrasound Cavitation, Obesity/Overweight

Induction of Fat Apoptosis: Mechanism of Action of Non-Invasive High-Intensity Electromagnetic Technology in a Porcine Model

Robert A. Weiss, MD, FAAD^{1*} and Jan Bernardy, MV²
¹Maryland Laser Skin, & Vein Institute, Hunt Valley, MD
²Veterinary Research Institute, Brno, Czech Republic

Objectives:

While controlled thermal changes in subcutaneous tissue have been used to trigger apoptosis of fat cells and have been proven clinically efficacious, another mechanism of electromagnetic stress suggests that fat apoptosis could be achieved by a non-thermal manner. This animal model study investigates the use of non-invasive high-intensity magnetic field device to induce apoptosis in fat cells.

Methods: Yorkshire pigs (N=2) received one treatment (30 minutes) in the abdominal area using a High-Intensity Focused Electromagnetic (HIFEM) device. Punch biopsy samples of fat tissue and blood samples were collected at the baseline, 1 and 8 hours after the treatment. Biopsy samples were sectioned and evaluated for the levels of apoptosis index (AI) by the TUNEL method. Statistical significance was examined using the tANOVA and Tukey test (α 5%). Biopsy samples were also assessed for molecular biomarkers. Blood samples were evaluated to determine changes related to fat and muscle metabolism.

Results: In treated pigs, a statistically significant increase in the apoptotic index (AI) (P=1.17E-4) was observed. Significant difference was found between AI at baseline (AI=18.75%) and 8-hours post-treatment (AI=35.95%). Serum levels of fat and muscle metabolism indicated trends (FFA -0.32 mmol L⁻¹, -28.1%; TG -0.24 mmol L⁻¹, -31.8%; Glycerol -5.68 mg L⁻¹, -34.8%; CK +67.1 µkat L⁻¹, +227.8%; LDH +4.9 µkat L⁻¹, +35.4%) suggesting that both adipose and muscle tissue were affected by HIFEM treatment. No adverse events were noted to skin and surrounding tissue.

Conclusions: Application of a high-intensity electromagnetic field in a porcine model results in adipocyte apoptosis. The analysis of serum levels suggests that HIFEM treatment influences fat and muscle metabolism. *Lasers Surg Med.* 51:47-53, 2019. © 2018 The Authors. *Lasers Surg Med and Medicine* Published by Wiley Periodicals, Inc

Key words: apoptosis; fat disruption; HIFEM; magnet technology; non-thermal

Original Article

Med Lasers 2016; 5(11):29-33
 pISSN 2287-8300 - eISSN 2288-0224

Impact of Contactless Apoptosis-Inducing Human Skin Surface and Subcutaneous Layer Histology: A Pilot Study

Bonched Goo¹
 Dae Suk Kim²

¹Naeum Dermatology and Aesthetic Clinic, Seoul, Korea
²Kangseon Gilim Dermatology Clinic, Seoul, Korea

Background and Objectives

Radiofrequency (RF) technology method to reduce subcutaneous to measure the changes of human subcutaneous fat layer temperature change in porcine adipocytes due to apoptosis-inducing RF device.

Materials and Methods

A single pig was treated with RF of 200 Watts. The skin was biopsied. A female human volunteer was treated with RF. The temperatures of the skin and subcutaneous fat were measured during the therapy.

Results

Skin biopsy specimens from the shape and size. Many of the adipocytes showed condensed chromatin and fragmented nuclei, indicating adipocyte apoptosis. In the human skin, the temperature of the skin maintained a temperature of 43°C.

Conclusion

The contactless selective RF device reduced subcutaneous fat in porcine skin. Apoptosis of skin and subcutaneous fat was maintained. Further clinical trials are needed to confirm safety.

Key words

Radiofrequency; Subcutaneous

Received June 15, 2016
 Accepted June 17, 2016

OPEN

Radiofrequency Heating and High-Intensity Focused Electromagnetic Treatment Delivered Simultaneously: The First Sham-Controlled Randomized Trial

Julene B. Samuels, M.D.
 Bruce Katz, M.D.
 Robert A. Weiss, M.D.

¹Louisville, Ky.; ²New York, N.Y.; and ³Hunt Valley, Md.

Background: Radiofrequency-based and high-intensity focused electromagnetic (HIFEM)-based devices have proved effective and safe for abdominal body shaping. Radiofrequency is known to reduce adipose tissue, whereas HIFEM treatment is effective for muscle definition. The authors investigated the efficacy of a novel device delivering synchronized radiofrequency and HIFEM treatment simultaneously for abdominal toning and fat reduction.

Methods: Seventy-two patients were enrolled and randomly divided into active (n = 48; age, 45.5 ± 13.0 years) and sham groups (n = 24; age, 44.6 ± 12.3 years). Both groups received three treatments on the abdomen once a week. The intensity in the active group was set to maximum tolerable level; in the sham group, the intensities were set to 5 percent. Ultrasound images were taken before treatment and at 1, 3, and 6 months after treatment to examine changes in subcutaneous fat and rectus abdominis muscle thickness. Digital photographs were taken, and satisfaction and therapy comfort were assessed.

Results: Ultrasound images of the active group at 1 month showed significant (p < 0.05) reduction in adipose tissue thickness by 20.5 percent (4.8 ± 2.6 mm), whereas rectus abdominis muscle thickness increased by 21.5 percent (2.0 ± 0.8 mm). At 3 months improved to 28.3 percent (7.6 ± 3.7 mm) and 24.2 percent (2.3 ± 0.9 mm), respectively. Improvements were maintained at 6 months after treatment in the active group, whereas the sham group showed no significant changes. Treatments were found to be comfortable. The active group showed higher satisfaction with outcomes.

Conclusion: Active treatment utilizing simultaneous application of radiofrequency and HIFEM therapy resulted in a significant increase in rectus abdominis thickness and subcutaneous fat reduction, exceeding previously published results for separate HIFEM and radiofrequency treatments. (*Plast Reconstr Surg.* 149:893e, 2022.)
CLINICAL QUESTION/LEVEL OF EVIDENCE: Therapeutic, II.



From private practice, the Juvu Skin and Laser Center, and Maryland Dermatology, Laser, Skin, & Vein Institute. Received for publication October 14, 2020; accepted August 04, 2021.

Presented at the 2020 Annual Meeting of the American Society for Dermatologic Surgery, held virtually, October 9 through 11, 2020.

This trial is registered under the name "PMS and Radiofrequency for Abdominal Toning and Reduction of Subcutaneous Fat," ClinicalTrials.gov identification no. NCT04587986 (<https://clinicaltrials.gov/ct2/show/NCT04587986>).

Copyright © 2022 The Author(s). Published by Wolters Kluwer Health, Inc. on behalf of the American Society of Plastic Surgeons. This is an open-access article distributed under the terms of the Creative Commons Attribution-NonCommercial-No Derivatives License 4.0 (CCBY-NC-ND), where it is permissible to download and share the work provided it is properly cited. The work cannot be changed in any way or used commercially without permission from the journal.

DOI: 10.1097/PRS.00000000000009030

Disclosure: Dr. Samuels, Dr. Katz, and Dr. Weiss are medical advisors for BTL Industries, Inc. No funding was received for this article.


Body Shaping: Build Waist and Abdominal Muscles, Say Goodbye to Sagging Fat

Massage




Multiple airbags apply and release pressure, simulating kneading and massage motions by human hands, making waist and abdominal muscles firmer and contours smoother.

RF (Radio Freque ncy)



High-frequency electromagnetic waves penetrate 3-5mm beneath the skin, generating heat to stimulate collagen regeneration and enhance muscle elasticity in the waist and abdomen.

Red Light




Red light penetrates the epidermis, accelerates waste metabolism, reduces pigmentation and wrinkles, and helps improve stretch marks and surgical scars.

Strengthen Your Core: Technology That Goes Beyond Fat Burning to Provide a “SPA” for Your Organs



PEMF

7.8Hz ultra-low-frequency PEMF cross-medium penetration technology delivers energy intervention up to 20cm deep into biological tissue. By regulating cell membrane potential and mitochondrial metabolic pathways, it significantly downregulates the expression of pro-inflammatory factors. It targets and alleviates various types of inflammation, such as enteritis, nephritis, prostatitis, gynecological inflammation, lumbar muscle strain, and lumbar disc herniation.



Vibration / Massage / Heat Therapy / RF Radio Frequency

Deep Activation of Microcirculation and Cellular Vitality: Efficiently clears metabolic waste and simultaneously alleviates inflammatory responses.

Intelligent Enhancement of Intestinal Rhythmic Movement: Optimizes nutrient conversion efficiency that helps restore a natural and comfortable state, and improves various intestinal issues.

One-Button Switch for Three Energy Fields: Unlock the Ultimate Solution for Core Health



Body Shaping Mode

7-in-1 Targeted Fat Burning

(PEMF, Heat Therapy, Ultrasound, Massage, Vibration, RF Radio Frequency, Red Light)

Function: Fat burning and body shaping, fitness support



Massage Mode

5D Deep Repair

(PEMF, Massage, Vibration, Heat Therapy, RF Radio Frequency)

Function: Strengthens the core and provides organ care



Heat Therapy Mode

2-Effect Gentle Care

(PEMF, Heat Therapy)

Function: Keeps the waist warm while supporting menstrual care

Two buttons for all operations

Shaken Massager Shaken Massager has simple function buttons, with a minimalist design of intuitive buttons, making it easy to switch modes and intensity.



- ① Long press the [⏻] key for 2 seconds to turn on the machine; when powered on, long press the [⏻] key for 2 seconds to turn off the machine.
- ② Short press the [⏻] key to switch usage levels.
- ③ Short press the [M] key to adjust modes (body shaping - massage - heat compress), with the three modes cycling in sequence; the product has a memory function and will default to the last mode before shutdown when powered on.
- ④ Long press the [M] key for 5 seconds to enter language switching, where the mode screen displays 8 languages for selection; short press the [M] key to choose a language, and short press the [⏻] key to confirm.

Wrist smart control, remote operation

Shaken Massager



Oly-Watch



Shaken Massager has simple function buttons and can be controlled via smartwatch touch, easily switching modes and intensity.

Built-in Bluetooth: quickly connect via Bluetooth and operate the device in real-time through the watch interface.

AI smart service, making every experience customized.

Under the AI smart algorithm, it automatically adjusts the working intensity and duration of technologies such as ultrasound, vibration, heat compress, and RF radio frequency.

It forms a total 20-minute operating plan intelligently combined from multiple sub-functions.



- 1、 Correctly wear it on the wrist.
- 2、 Connect to Shaken Massager.



Enter AI mode. Collect user body data: heart rate, blood oxygen, sleep quality



Using an exclusive AI algorithm, it intelligently generates a personalized working program that best fits the user's current body state




Experience the personalized fat-shaking service provided by the product, making it easy, efficient, and worry-free.

A thousand people, a thousand bodies—AI tailors it for you.

Wear Freely, Burn Fat with Elegance


Dual-Screen Control at Your Fingertips

Real-time display of operational status: the left screen shows gear and intensity, while the right screen displays mode and countdown timer.




Seamless Fit to Your Waistline

Featuring a joint with **70°-120°** degrees of freedom for adjustment, it accommodates waistlines from 60 to 150 cm and is suitable for individuals weighing from 40 to 130 kg, ensuring it does not slip when bending or turning.



Technology That Balances Strength and Comfort


Constructed from soft silicone material and precision cast 304 stainless steel powder, it exhibits stable properties, good thermal conductivity, wear and corrosion resistance, is hypoallergenic, and can be easily cleaned with an alcohol wipe.





Olywatch

sleeve



Unboxing complete contents

**The belts on both sides of the equipment
cannot be disassembled
Wearing sleeve covers makes cleaning convenient**

Waist Extension Strap



Correct installation of the adhesive Strap



Tail end folded outward ✓

Incorrect installation of the adhesive Strap



Tail end folded inward ✗

Suitable for



Sedentary individuals (to relieve lower back muscle strain)



Fitness and body-shaping enthusiasts



Individuals with abdominal obesity



Postpartum women (to improve the appearance of stretch marks)



Individuals with intestinal disorders



Those experiencing dysmenorrhea, irregular menstruation, or gynecological inflammation




Individuals with prostatitis, urinary urgency, or incomplete urination



Those with nephritis, cystitis, or similar conditions

Tech-Powered Fat-Burning, Easy and Comfortable




For professional use



For recreational use



For athletic activities



For business trips



Shaken Massager

PEMF Smart Fat-Burning Device

Shake off fat and shine with confidence every day!